

Canova

ERGO design
for greater user comfort



Canova is an embedding module designed to meet operator needs during histological specimen embedding.



Ergonomics: innovative philosophy ERGO-TECH assures to the operator the best comfort and ergonomics with rounded edges and comfortable soft coverings.



Safety: the integrated suction system prevents operator's exposure to fumes and paraffin vapors. The FMC option controls and manages the periodic maintenances to guarantee always quality standard.



Time saving: patented system for fast paraffin melting (about 3 hours). Paraffin collecting drawer with disposable container made of plasticized paper.

More new features



Innovative working plane

Equipped with a heated grid plate to clean cassettes from paraffin in excess.



Canova Maintenance Kit

A guarantee for instrument efficiency and operator safety.



Heated tamper and Canova Safe Forceps

Manageable by software with hour countdown control system.

Canova cooling plates

Cooling device arranged in two plates for a better space management. Capacity up to 110 standard blocks.



For further information, please reply directly to the present e-mail or contact:

Mr. Marco Badolato

☎ (+39) 0363-986411

✉ (+39) 0363-948000

✉ marco.badolato@diapath.com